

AD-A193 264


MUTAGENIC POTENTIAL OF NITROGUANIDINE IN THE AMES  
SALMONELLA/MAMMALIAN MI. (U) LETTERMAN ARMY INST OF  
RESEARCH PRESIDIO OF SAN FRANCISCO CA.  
S E SEBASTIN ET AL. MAR 88 LAIR-260


1/1

UNCLASSIFIED

F/G 6/5


NL





MICROCOPY RESOLUTION TEST CHART  
 NBS 1963-A

**AD-A193 264**



UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE

## REPORT DOCUMENTATION PAGE

Form Approved  
OMB No. 0704-0188

1a. REPORT SECURITY CLASSIFICATION UNCLASSIFIED		1b. RESTRICTIVE MARKINGS	
2a. SECURITY CLASSIFICATION AUTHORITY		3. DISTRIBUTION / AVAILABILITY OF REPORT Approved for public release; distribution is unlimited.	
2b. DECLASSIFICATION / DOWNGRADING SCHEDULE			
4. PERFORMING ORGANIZATION REPORT NUMBER(S) IAIR Institute Report No. 260		5. MONITORING ORGANIZATION REPORT NUMBER(S)	
6a. NAME OF PERFORMING ORGANIZATION Genetic Toxicology Branch Division of Toxicology	6b. OFFICE SYMBOL (If applicable) SGRD-UL-TO-G	7a. NAME OF MONITORING ORGANIZATION US Army Biomedical Research and Development Laboratory	
6c. ADDRESS (City, State, and ZIP Code) Letterman Army Institute of Research Presidio of San Francisco, CA 94129-6800		7b. ADDRESS (City, State, and ZIP Code) Ft. Detrick Frederick, MD 21701-5010	
8a. NAME OF FUNDING / SPONSORING ORGANIZATION US Army Medical Research and Development Command	8b. OFFICE SYMBOL (If applicable) SGRD-2A	9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER	
8c. ADDRESS (City, State, and ZIP Code) Fort Detrick Frederick, MD 21701-5012		10. SOURCE OF FUNDING NUMBERS	
		PROGRAM ELEMENT NO 62720A	PROJECT NO 835
		TASK NO AB	WORK UNIT ACCESSION NO. DA 303913
11. TITLE (Include Security Classification) Mutagenic potential of nitroguanidine in the Ames <u>Salmonella</u> /mammalian microsome mutagenicity test			
12. PERSONAL AUTHOR(S) Suzanne E. Sebastian and Don W. Korte, Jr			
13a. TYPE OF REPORT	13b. TIME COVERED FROM _____ TO _____	14. DATE OF REPORT (Year, Month, Day) March 1988	15. PAGE COUNT 17
16. SUPPLEMENTARY NOTATION			
17. COSATI CODES		18. SUBJECT TERMS (Continue on reverse if necessary and identify by block number)	
FIELD	GROUP	SUB-GROUP	
19. ABSTRACT (Continue on reverse if necessary and identify by block number) The mutagenic potential of NITROGUANIDINE was assessed by using the Ames <u>Salmonella</u> /Mammalian Microsome Mutagenicity Test. Tester strains TA97, TA98, TA100, TA102, TA1535, TA1537, and TA1538 were exposed to doses ranging from 2.8 mg/plate to 0.0875 mg/plate. The test compound was not mutagenic under conditions of this test.			
20. DISTRIBUTION / AVAILABILITY OF ABSTRACT <input checked="" type="checkbox"/> UNCLASSIFIED/UNLIMITED <input type="checkbox"/> SAME AS RPT <input type="checkbox"/> DTIC USERS		21. ABSTRACT SECURITY CLASSIFICATION unclassified	
22a. NAME OF RESPONSIBLE INDIVIDUAL Edwin S. Beatrice, COL, MC		22b. TELEPHONE (Include Area Code) (415) 561-3600	22c. OFFICE SYMBOL SGRD-UL-2

# ABSTRACT

The mutagenic potential of NITROGUANIDINE was assessed by using the Ames Salmonella/Mammalian Microsome Mutagenicity Test. Tester strains TA97, TA98, TA100, TA102, TA1535, TA1537, and TA1538 were exposed to doses ranging from 2.8 mg/plate to 0.0875 mg/plate. The test compound was not mutagenic under conditions of this test.

Key Words: Mutagenicity, Genetic Toxicology, Ames Test, NITROGUANIDINE.



Accession For	
NTIS GRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A-1	

## PREFACE

TYPE REPORT: Ames Test GLP Study Report

TESTING FACILITY: US Army Medical Research and Development  
Command  
Letterman Army Institute of Research  
Presidio of San Francisco, CA 94129-6800

SPONSOR: US Army Medical Research and Development Command  
US Army Biomedical Research and Development  
Laboratory  
Fort Detrick, Frederick, MD, 21701-5010  
Project Officer: Gunda Reddy, PhD

PROJECT/WORK UNIT/APC: 3E162720A835/180/TLBO

GLP STUDY NUMBER: 86008

STUDY DIRECTOR: MAJ Don W. Korte Jr, PhD, MSC

PRINCIPAL INVESTIGATOR: Suzanne E. Sebastian, BA, SP4, USA

REPORT AND DATA MANAGEMENT: A copy of the final report,  
study protocol, retired SOP's,  
stability and purity data on the  
test compound, and an aliquot of the  
test compound will be retained in  
the LAIR Archives.

TEST SUBSTANCE: NITROGUANIDINE

INCLUSIVE STUDY DATES: 28-30 October 1986

OBJECTIVE: The objective of this study was to determine the  
mutagenic potential of NITROGUANIDINE (LAIR Code TP 36A) by  
using the Ames Salmonella/Mammalian Microsome Mutagenicity  
Test.


#### ACKNOWLEDGMENTS


MAJ John W. Harbell, PhD, MSC; SGT Lillie D. Witcher, BS; and SGT Gayle A. Orner, BS, provided research assistance.



SIGNATURES OF PRINCIPAL SCIENTISTS INVOLVED IN THE STUDY

We, the undersigned, declare that GLP Study 86008 was performed under our supervision, according to the procedures described herein, and that this report is an accurate record of the results obtained.

  
DON W. KORTE Jr, PhD / Date  
MAJ, MSC  
Study Director

  
SUZANNE E. SEBASTIAN, BA / DATE  
SP4, USA  
Principal Investigator



DEPARTMENT OF THE ARMY

LETTERMAN ARMY INSTITUTE OF RESEARCH  
PRESIDIO OF SAN FRANCISCO, CALIFORNIA 94129-6800

REPLY TO  
ATTENTION OF:

SGRD-ULZ-QA (70-1n)

23 February 1988

MEMORANDUM FOR RECORD

SUBJECT: Report of GLP Compliance for 86008, Tox Series 107

1. I hereby certify that the protocol was reviewed on 29 October 1986.
2. The report and raw data for this study were audited on 24 November 1987.

*Carolyn M. Lewis*  
CAROLYN M. LEWIS  
C, Quality Assurance

## TABLE OF CONTENTS

Abstract .....	i
Preface .....	ii
Acknowledgments.....	iii
Signatures of Principal Scientists.....	iv
Report of the Quality Assurance Unit.....	v
Table of Contents.....	vi
BODY OF THE REPORT	
INTRODUCTION	
Objective of the Study .....	1
MATERIALS AND METHODS	
Test Compound .....	2
Test Solvent .....	2
Chemical Preparation .....	2
Test Strains .....	2
Mammalian Microsome System .....	3
Test Format .....	3
Data Interpretation .....	4
Deviations from the Protocol/SOP .....	5
Storage of the Raw Data and Final Report .....	5
RESULTS.....	5
DISCUSSION.....	9
CONCLUSION.....	9
REFERENCES.....	10

APPENDIX A.....	11
APPENDIX B.....	14
OFFICIAL DISTRIBUTION LIST.....	18

**Mutagenic Potential of NITROGUANIDINE - Sebastian and Korte**

NITROGUANIDINE, a primary component of US Army triple-base propellants, is now produced in a Government-owned contractor-operated ammunition plant. The US Army Biomedical Research and Development Laboratory (USABRDL), as part of its charge to evaluate the environmental and health hazards of propellants generated by US Army munitions manufacturing facilities, conducted a review of the nitroguanidine data base and identified significant gaps in the toxicity data (1). The Division of Toxicology, LAIR, was tasked by USABRDL to develop a genetic and mammalian toxicity profile for nitroguanidine and related intermediates/by-products of its manufacture or environmental degradation products.

The Ames Salmonella/Mammalian Microsome Mutagenicity Test is a short-term screening test that utilizes histidine auxotrophic mutant strains of Salmonella typhimurium to detect compounds that are potentially mutagenic in mammals. A mammalian microsomal enzyme system is incorporated in the test to increase sensitivity by simulating in vivo metabolic activation of the test compound. The Ames test is an inexpensive yet highly predictive and reliable test for detecting mutagenic activity and thus carcinogenic potential (2).

This evaluation of NITROGUANIDINE utilizes a revision of the Ames Salmonella/Mammalian Microsome Mutagenicity Test (3). Two new tester strains, a frame-shift strain (TA97) and a strain carrying an ochre mutation on a multicopy plasmid (TA102), are added to the standard tester set.

Sebastian and Korte--2

### Objective of the Study

The objective of this study was to determine the mutagenic potential of NITROGUANIDINE (LAIR Code TP 36A) by using the revised Ames Salmonella/Mammalian Microsome Mutagenicity Test.

### **MATERIALS AND METHODS**

#### Test Compound

Chemical name: NITROGUANIDINE

Code number: LAIR Code No. TP 36A

Physical state: White crystalline solid

Source: Sunflower Army Ammunition Plant  
De Soto, KS

Storage: NITROGUANIDINE was received from Sunflower Army Ammunition Plant, De Soto, KS, and assigned the LAIR Code number TP 36A. The test compound was stored at room temperature (21°C) until used.

Chemical Properties/Analysis: Data provided by Sunflower Army Ammunition Plant characterizing the chemical composition and purity of the test material are presented in Appendix A with confirmatory analysis of the test material performed by the Division of Toxicology, LAIR (Presidio of San Francisco, CA).

#### Test Solvent

The positive control chemicals were dissolved in grade I dimethyl sulfoxide (lot 113F-0450) obtained from Sigma Chemical Co. (St. Louis, MO). The test chemical was dissolved in the same lot of DMSO. Reagent grade water used in this assay is first passed through a Technic Series 300 Reverse Osmosis Unit (Seattle, WA), then through a Corning MP-1 Mega Pure System glass distillation unit (Corning Glass Works, Corning, NY) (4).

#### Chemical Preparation

On the day of dosing, the compound was dissolved directly into DMSO at a concentration of 28 mg/ml for the highest dose. Aliquots of this solution were used to prepare the serial dilutions.

### Test Strains

Salmonella strains TA97, TA98, TA100, TA102, TA1535, TA1537, and TA1538 obtained from the laboratory of Dr. Bruce Ames, University of California, Berkeley, were used. These strains were maintained in our laboratory at -80°C. Quality control tests were run concurrently with the test substance to establish the validity of their special features and to determine the spontaneous reversion rate. Descriptions of the strains, their genetic markers, and the methods for strain validation are given in the LAIR SOP, OP-STX-1 (5).

### Mammalian Microsome System

The S-9 (batch #R-315) was purchased from Microbiological Associates Inc. (Bethesda, MD). The optimal titer of this S-9, as determined by Microbiological Associates Inc., was 0.75 mg protein/plate.

### Test Format

NITROGUANIDINE was evaluated for mutagenic potential according to the revised Ames method (3). A detailed description of the methodology is given in LAIR SOP, OP-STX-1 (5). The preincubation modification was chosen to enhance the sensitivity of the assay by exposing the bacteria to higher concentrations of test compound (and the activation products, when present) than would be possible in the standard plate incorporation assay. The bacteria were preincubated in the presence of the compound, both with and without metabolic activation, for 20 minutes on a shaker incubator at 37°C. A single preincubation tube was prepared for each top agar triplicate. Each preincubation tube contained 10 ml of a mixture which consisted of 1 ml of bacteria (16 hour culture), 1 ml of test compound (28 mg/ml or a serial dilution), 2 ml of S-9 if required, and the remaining volume nutrient broth. The highest dose (2.8 mg/ml) in the preincubation mixture approached the practical limits for nitroguanidine solubility in aqueous media (1). The top agar tubes were prepared by adding 0.7 ml from the preincubation tube to 2 ml of top agar. After mixing, the top agar was then overlaid on minimal glucose agar plates. These plates contained 2% glucose and Vogel Bonner "E" concentrate (6). Plates were incubated upside down in the dark at 37°C for 72 hours (Maron, 1985, personal communication). Plates were prepared in triplicate and the individual revertant counts were recorded. The average number of revertants at each dose level was compared to the average number of spontaneous revertants (negative control). The spontaneous reversion rate (with and without S-9) was

monitored by averaging the counts from two determinations run simultaneously with the test compound. The spontaneous reversion rate was determined by inoculating one set of plates before and one set after the test compound plates so that any change in spontaneous reversion rate during the dosing procedure would be detected. This spontaneous reversion rate was also compared with historical values for this laboratory and those cited in Maron and Ames (3). Sterility and strain verification controls were run concurrently. All reagents, test compounds, and media were checked for sterility by plating samples of each on MGA media and incubating them at 37°C with the test plates. The Salmonella strains were verified by a standard battery of tests. The integrity of the different Salmonella strains used in the assay was verified by the following standard tests:

- Lack of growth (inhibition) in the presence of crystal violet which indicated that the prerequisite alteration of the lipopolysaccharide layer of the cell wall was present.

- Growth in the presence of ampicillin-impregnated disks which indicated the presence of an ampicillin-resistant R Factor in all strains except TA1535, TA1537, and TA1538.

- Lack of growth (inhibition) following exposure to ultraviolet light which indicated the absence of the DNA excision-repair mechanism (for all strains except TA102).

Six known mutagens were tested as positive controls to confirm the responsiveness of the strains to the mutation process. Each strain must be tested with at least one positive control but may be tested with several. These compounds, benzo[a]pyrene (lot 18C-0378), 2-aminofluorene (lot 021547), 2-aminoanthracene (lot 020797), mitomycin-C (lot 015F-0655), 4-nitroquinoline-n-oxide (lot 89C-0710) and N-methyl-N'-nitro-N-nitrosoguanidine (lot 127C-0342), were obtained from Sigma Chemical Co. (St. Louis, MO). The test compound and mutagens were handled during this study in accordance with the standards published in NIH Guidelines for the Laboratory Use of Chemical Carcinogens (DHHS Publication No. (NIH) 81-2385, May 1981).

#### Data Interpretation

According to Brusick (7), a compound is considered mutagenic if a positive dose response (correlated dose



response) over three dose concentrations is achieved with at least the highest dose yielding a revertant colony count greater than or equal to twice the spontaneous colony count for the tester strains TA98 and TA100, or three times the spontaneous colony count for strains TA1535, TA1537, and TA1538 (3,5). A strong correlated dose response in strain TA100 without a doubling of the individual colony count may also be considered positive.

Maron and Ames (3) consider a compound mutagenic in tester strains TA97 and TA102 if a correlated dose response over three concentrations is achieved with the highest dose yielding a revertant colony count greater than or equal to twice the spontaneous colony count.

#### Deviations from the Protocol/SOP

The preincubation modification of the Ames Assay was chosen to enhance sensitivity by exposing bacteria to a higher concentration of compound, for a longer period of time. Volumes for the preincubation mixture were different from those specified in the SOP because of the limited solubility of nitroguanidine. This deviation has no impact on the validity of the study.

#### Storage of the Raw Data and Final Report

A copy of the final report, study protocols, raw data, SOPs, and an aliquot of the test compound will be retained in the LAIR archives.

### **RESULTS**

Normal results were obtained for all sterility and strain verification tests during the Ames Test performed on 21-24 May, 1986 (Table 1). NITROGUANIDINE did not induce any appreciable increase in the revertant colony counts relative to those of the negative control cultures (Table 2).

A tabular presentation of raw data is included in Appendix B.

TABLE 1

STRAIN VERIFICATION AND STERILITY TESTING FOR  
THE MUTAGENICITY DETERMINATION  
ON NITROGUANIDINE (TP 36A)

GLP STUDY NUMBER 86008

STRAIN VERIFICATION

OBSERVATIONS\*

STRAIN	HISTIDINE REQUIREMENT	AMPICILLIN RESISTANCE	UV REPAIR	CRYSTAL VIOLET	STERILITY CONTROL
TA97	NG	G	NG	NG	NG
TA98	NG	G	NG	NG	NG
TA100	NG	G	NG	NG	NG
TA102	NG	G	G	NG	NG
TA1535	NG	NG	NG	NG	NG
TA1537	NG	NG	NG	NG	NG
TA1538	NG	NG	NG	NG	NG

STERILITY CONTROL FOR MUTAGENICITY DETERMINATION

<u>MATERIAL TESTED</u>	<u>OBSERVATION*</u>
MINIMAL GLUCOSE AGAR PLATES	NG
TOP AGAR	NG
DILUENT WATER	NG
NUTRIENT BROTH	NG
TEST COMPOUND (HIGHEST DOSE)	NG
S-9	NG

---

\*G = Growth, NG = No Growth

TABLE 2

## NITROGUANIDINE

## REVERTANTS/PLATE

MEAN  $\pm$  1SD

COMPOUND*	DOSE/PLATE	TA97	TA98	TA100	TA102
<b>WITHOUT S-9</b>					
NEG CONTROL	0.0 mg	88 $\pm$ 16.1	14 $\pm$ 3.3	111 $\pm$ 10.2	173 $\pm$ 13.2
MITO C	0.5 $\mu$ g				1327 $\pm$ 22.5
MNNG	2.0 $\mu$ g			1460 $\pm$ 149.1	
MNNG	20.0 $\mu$ g				
NQNO	2.0 $\mu$ g	252 $\pm$ 92.2			
TP 36A	2.8 mg	66 $\pm$ 10.6	19 $\pm$ 0.6	88 $\pm$ 22.5	119 $\pm$ 1.7
TP 36A	1.4 mg	96 $\pm$ 14.6	13 $\pm$ 3.8	107 $\pm$ 2.9	198 $\pm$ 2.0
TP 36A	0.7 mg	83 $\pm$ 15.3	12 $\pm$ 4.6	106 $\pm$ 11.0	166 $\pm$ 24.4
TP 36A	0.35 mg	102 $\pm$ 27.7	12 $\pm$ 3.2	85 $\pm$ 15.5	162 $\pm$ 15.0
TP 36A	0.175 mg	96 $\pm$ 7.6	20 $\pm$ 2.0	96 $\pm$ --	146 $\pm$ 35.6
TP 36A	0.0875 mg	93 $\pm$ 6.1	16 $\pm$ 5.5	99 $\pm$ 13.1	172 $\pm$ 12.7
<b>WITH S-9</b>					
NEG CONTROL	0.0 mg	110 $\pm$ 16.6	35 $\pm$ 10.8	125 $\pm$ 21.7	282 $\pm$ 5.4
2-AA	2.0 $\mu$ g		298 $\pm$ 73.4	563 $\pm$ 86.5	
2-AF	2.0 $\mu$ g	572 $\pm$ 136.7	1683 $\pm$ 42.5	706 $\pm$ 86.5	
BP	2.0 $\mu$ g	251 $\pm$ 7.5	125 $\pm$ 7.5	366 $\pm$ 35.8	
TP 36A	2.8 mg	144 $\pm$ 38.3	38 $\pm$ 8.5	110 $\pm$ 36.2	272 $\pm$ 27.2
TP 36A	1.4 mg	107 $\pm$ 9.5	25 $\pm$ 2.1	122 $\pm$ 26.6	282 $\pm$ 20.4
TP 36A	0.7 mg	121 $\pm$ 7.6	30 $\pm$ 1.7	99 $\pm$ 6.4	308 $\pm$ 13.9
TP 36A	0.35 mg	101 $\pm$ 9.0	26 $\pm$ 1.0	127 $\pm$ 22.7	280 $\pm$ 10.2
TP 36A	0.175 mg	115 $\pm$ 4.6	27 $\pm$ 0.6	105 $\pm$ --	267 $\pm$ 59.5
TP 36A	0.0875 mg	107 $\pm$ 6.6	33 $\pm$ 4.5	129 $\pm$ 9.7	298 $\pm$ 15.1

\*MITO-C=mitomycin-C, MNNG=N-methyl-N'-nitro-N-nitrosoguanidine, NQNO=4-nitroquinoline-n-oxide, AA=2-aminoanthracene, AF=2-aminofluorene, BP=benzo[a]pyrene.

TABLE 2, (continued)

## NITROGUANIDINE

## REVERTANTS/PLATE

MEAN  $\pm$  1SD

COMPOUND*	DOSE/PLATE	TA1535	TA1537	TA1538
<b>WITHOUT S-2</b>				
NEG CONTROL	0.0 mg	11 $\pm$ 3.0	10 $\pm$ 3.2	18 $\pm$ 6.4
MNNG	20.0 $\mu$ g	2224 $\pm$ 378.2		
TP 36A	2.8 mg	21 $\pm$ 2.5	19 $\pm$ 4.6	19 $\pm$ 4.0
TP 36A	1.4 mg	13 $\pm$ 4.7	7 $\pm$ 0.6	10 $\pm$ 4.6
TP 36A	0.7 mg	10 $\pm$ 2.3	12 $\pm$ 5.6	9 $\pm$ 3.6
TP 36A	0.35 mg	10 $\pm$ 1.2	10 $\pm$ 0.0	16 $\pm$ 5.9
TP 36A	0.175 mg	14 $\pm$ 4.0	6 $\pm$ 2.0	13 $\pm$ 0.6
TP 36A	0.0875 mg	16 $\pm$ 6.4	5 $\pm$ 0.6	17 $\pm$ 1.0
<b>WITH S-2</b>				
NEG CONTROL	0.0 mg	11 $\pm$ 1.7	12 $\pm$ 4.1	27 $\pm$ 2.4
2-AA	2.0 $\mu$ g		318 $\pm$ 55.2	572 $\pm$ 231.1
2-AF	2.0 $\mu$ g			1168 $\pm$ 173.9
BP	2.0 $\mu$ g		48 $\pm$ 5.9	74 $\pm$ 39.0
TP 36A	2.8 mg	9 $\pm$ 3.6	10 $\pm$ 3.6	25 $\pm$ 9.6
TP 36A	1.4 mg	18 $\pm$ 13.0	10 $\pm$ 4.2	34 $\pm$ 1.7
TP 36A	0.7 mg	12 $\pm$ 4.9	9 $\pm$ 2.5	25 $\pm$ 4.9
TP 36A	0.35 mg	14 $\pm$ 4.5	11 $\pm$ 4.2	14 $\pm$ 6.5
TP 36A	0.175 mg	15 $\pm$ 3.0	12 $\pm$ 4.6	24 $\pm$ 7.0
TP 36A	0.0875 mg	16 $\pm$ 11.8	12 $\pm$ 2.5	17 $\pm$ 1.0

\*MITO-C=mitomycin-C, MNNG=N-methyl-N'-nitro-N-nitrosoguanidine, NQNO=4-nitroquinoline-n-oxide, AA=2-aminoanthracene, AF=2-aminofluorene, BP=benzo[a]pyrene.

## DISCUSSION

Certain test criteria must be satisfied before an Ames test can be considered a valid assessment of a compound's mutagenic potential. First, the special features of the Ames strains must be verified. These features include demonstration of ampicillin resistance, alterations in the lipopolysaccharide layer, and deficiency in DNA excision-repair (except TA102). Second, the Salmonella strains must be susceptible to mutation by known mutagens. Third, the optimal concentration of the test compound must be determined by treating TA100 with a broad range of doses and observing the potential toxic effects on formation of macrocolonies and microcolonies. If these tests are performed and expected data are obtained, then the results of an Ames test can be considered valid.

After validation of bacterial strains and selection of optimal sublethal doses, NITROGUANIDINE was evaluated in the Ames preincubation test. Criteria for a positive response include both a correlated dose response over three dose concentrations, and a revertant colony count at least two times for TA97, TA98, TA100, and TA102 (2,7) or three times for TA1535, TA1537, and TA1538 (3,5) the spontaneous revertant colony count. NITROGUANIDINE did not induce the requisite dose-response relationship or the increase in revertant colony counts necessary for a positive response. Thus, the results of this test indicate that NITROGUANIDINE is not mutagenic when evaluated in the Ames test.

## CONCLUSION

NITROGUANIDINE was evaluated for mutagenic potential in the Ames test, both in the presence and absence of metabolic activation, and did not induce a positive mutagenic response under conditions of this study.

## REFERENCES

1. Kenyon KF. A data base assessment of environmental fate aspects of nitroguanidine. Frederick, Maryland: US Army Bioengineering Research and Development Laboratory, 1982, DTIC No. ADA 125591.
2. Ames BN, McCann J, Yamasaki E. Methods for detection of carcinogens and mutagens with Salmonella/Mammalian Microsome Mutagenicity Test. *Mutat Res* 1975;31:347-364.
3. Maron DM, Ames BN. Revised methods for the Salmonella Mutagenicity Test. *Mutat Res* 1983;113:173-215.
4. Operation of the Technic Model 301 Reverse Osmosis Pre-Treatment Water System and the Corning Model MP-1 Glass Still. LAIR Standard Operating Procedure OP-STX-94, Presidio of San Francisco, California: Letterman Army Institute of Research, 29 July 1985.
5. Ames Salmonella/Mammalian Microsome Mutagenesis Test. LAIR Standard Operating Procedure OP-STX-1, Presidio of San Francisco, California: Letterman Army Institute of Research, 29 August 1986.
6. Vogel HJ, Bonner DM. Acetylornithinase of E. coli: Partial purification and some properties. *J Biol Chem* 1956;218:97-106.
7. Brusick D. Genetic toxicology. In: Hayes AW, ed. Principles and methods of toxicology. New York: Raven Press, 1982:223-272.

APPENDIX A


CHEMICAL DATA

Chemical name: Nitroguanidine (NGu)

Other listed names: Guanidine, Nitro; alpha-Nitroguanidine;  
beta-Nitroguanidine

LAIR Code: TP 36A

Structural formula:



Molecular formula: CH4N4O2

Molecular weight: 104.1

pH range of dosing suspensions: 6.7 - 7.4(1)

Physical state: White Powder

Melting point: 232° C(2)

Source: Hercules Aerospace Division  
Sunflower Ammunition Plant  
DeSoto, Kansas

Lot No. SOW84K010-A-001

Analytical data/purity:

The major peaks in the infrared spectrum of the compound were observed at 3450, 3396, 3342, 3278, 3201, 1666, 1634, 1525, 1404, 1314, 1151, 1045, 782 cm<sup>-1</sup>. (3) The spectrum obtained for the test compound in our lab was identical to the Sadtler standard spectrum for nitroguanidine. (4) HPLC showed only one peak (retention time 4.9 min). (5) The conditions employed were as follows: column, Brownlee RP-18 (4.6 x 250 mm); solvent, 10% methanol-90% water; flow rate, 0.7 ml/min; oven temperature, 50°C; monitoring wavelength, 265 nm.

1. Wheeler CR. Nitrocellulose-Nitroguanidine Projects  
Laboratory Notebook #85-12-022, p 26. Letterman Army  
Institute of Research, Presidio of San Francisco, CA.

Sebastian and Korte--12

2. Fedoroff BT, Sheffield OE. Encyclopedia of explosives and related items. Vol V. Dover, New Jersey: Picatinny Arsenal 1975: G154.
3. Wheeler, CR. Nitrocellulose-Nitroguanidine Projects. Laboratory Notebook #85-12-022, p. 22-23. Letterman Army Institute of Research, Presidio of San Francisco, CA.
4. Sadtler Research Laboratory, Inc. Sadtler standard spectra. Philadelphia: The Sadtler Research Laboratory, Inc., 1962: Infra-red spectrogram #21421.
5. Wheeler, CR. Nitrocellulose-Nitroguanidine Projects Laboratory Notebook #85-12-022, pp. 24-25. Letterman Army Institute of Research, Presidio of San Francisco, CA.



APPENDIX B-1  
RAW DATA TABLE  
NEGATIVE CONTROL  
NITROGUANIDINE (TP 36A)

COMPOUND	DOSE	TA97	TA98	TA100	TA102	TA1535	TA1537	TA1538
<b><u>WITHOUT S-9</u></b>								
NEG CONTROL (START RUN)	0.0 mg/plate	101 79 61	14 12 10	105 107 97	178 193 174	16 8 11	11 14 13	21 24 26
NEG CONTROL (END RUN)	0.0 mg/plate	100 101 85	12 16 19	115 126 117	158 159 178	12 10 8	8 6 8	13 10 15
<b><u>WITH S-9</u></b>								
NEG CONTROL (START RUN)	0.0 mg/plate	104 92 100	45 36 46	134 116 99	283 275 280	13 9 13	10 18 6	26 29 27
NEG CONTROL (END RUN)	0.0 mg/plate	140 113 113	36 29 17	149 * *	288 * *	12 10 10	15 12 11	28 22 27

\*plate contaminated

APPENDIX B-2  
RAW DATA TABLE  
POSITIVE CONTROLS

COMPOUND*	DOSE	TA97	TA98	TA100	TA102	TA1535	TA1537	TA1538
2-AA	2.0 µg/plate		377	565			376	566
			285	476			266	806
			232	649			313	344
2-AF	2.0 µg/plate	730	1683	682				1235
		501	1641	634				1299
		486	1726	802				971
BP	2.0 µg/plate	244	118	335			46	43
		259	133	405			55	62
		251	125	357			44	118
MITO C	0.5 µg/plate				1341			
					1301			
					1339			
MNNG	2.0 µg/plate			1601				
				1304				
				1475				
MNNG	20.0 µg/plate					1801		
						2340		
						2530		
NQNO	2.0 µg/plate	289						
		147						
		320						

\*MITO-C=mitomycin-C, MNNG=N-methyl-N'-nitro-N-nitrosoguanidine, NQNO=4-nitroquinoline-n-oxide, AA=2-aminoanthracene, AF=2-aminofluorene, BP=benzo[a]pyrene.

## APPENDIX B-3

## RAW DATA TABLE

## NITROGUANIDINE (TP 36A)

NITROUT S-2

COMPOUND	DOSE	TA97	TA98	TA100	TA102	TA1535	TA1537	TA1538
TP 36A	2.8 mg/p	68	19	80	118	21	24	23
		55	18	70	121	24	15	20
		76	19	113	118	19	18	15
TP 36A	1.4 mg/p	84	15	109	198	8	7	15
		112	9	109	196	15	6	9
		91	16	104	200	17	7	6
TP 36A	0.7 mg/p	89	15	93	139	7	17	10
		95	7	113	174	11	13	5
		66	15	111	186	11	6	12
TP 36A	0.35 mg/p	133	10	100	177	11	10	18
		95	16	86	147	9	10	9
		79	11	69	162	11	10	20
TP 36A	0.175 mg/p	87	22	96	180	16	8	13
		101	20	*	150	9	4	14
		99	18	*	109	16	6	13
TP 36A	0.0875 mg/p	100	16	111	186	9	5	16
		90	22	85	161	21	5	17
		89	11	100	170	19	4	18

\*plate contaminated

## APPENDIX B-4

## RAW DATA TABLE

## NITROGUANIDINE (TP 36A)

WITH S-9

COMPOUND	DOSE	TA97	TA98	TA100	TA102	TA1535	TA1537	TA1538
TP 36A	2.8 mg/p	103 179 149	37 30 47	92 152 87	303 254 258	8 13 6	11 13 6	27 15 34
TP 36A	1.4 mg/p	110 96 114	27 24 23	106 108 153	259 297 291	10 11 33	9 7 15	36 33 33
TP 36A	0.7 mg/p	114 129 119	31 28 31	94 * 103	300 324 300	8 15 *	9 11 6	23 22 31
TP 36A	0.35 mg/p	95 96 111	27 25 26	119 153 110	287 268 284	10 19 14	10 8 16	14 7 20
TP 36A	0.175 mg/p	120 114 111	27 28 27	105 * *	200 313 289	12 18 16	17 9 9	19 21 32
TP 36A	0.0875 mg/p	100 108 113	29 33 38	118 131 137	296 284 314	9 10 30	15 12 10	18 17 16

\*plate lost

## Distribution List

Commander  
US Army Biomedical Research and  
Development Laboratory (27)  
ATTN: SGRD-UBZ-C  
Fort Detrick, Frederick, MD 21701-5010

Defense Technical Information Center  
(DTIC) (2)  
ATTN: DTIC-DLA  
Cameron Station  
Alexandria, VA 22304-6145

US Army Medical Research and  
Development Command (2)  
ATTN: SGRD-RMI-S  
Fort Detrick, Frederick, MD 21701-5012

Commandant  
Academy of Health Sciences, US Army  
ATTN: AHS-CDM  
Fort Sam Houston, TX 78234

Chief  
USAEHA Regional Division, West  
Fitzsimmons AMC  
Aurora, CO 80045

Chief  
USAEHA Regional Division, North  
Fort George G. Meade, MD 20755

Chief  
USAEHA Regional Division, South  
Bldg. 180  
Fort McPherson, GA 30330

Commander  
USA Health Services Command  
ATTN: HSPA-P  
Fort Sam Houston, TX 78234

Commandant  
Academy of Health Sciences  
United States Army  
ATTN: Chief, Environmental  
Quality Branch  
Preventive Medicine Division  
(HSHA-IPM)  
Fort Sam Houston, TX 78234

Commander US Army Materiel  
Command  
ATTN: AMSCG  
5001 Eisenhower Avenue  
Alexandria, VA 22333

Commander  
US Army Environmental Hygiene  
Agency  
ATTN: Librarian, HSDH-AD-L  
Aberdeen Proving Ground, MD 21010

Dean  
School of Medicine  
Uniformed Services University of the  
Health Sciences  
4301 Jones Bridge Road  
Bethesda, MD 20014

Commander  
US Army Materiel Command  
ATTN: AMCEN-A  
5001 Eisenhower Avenue  
Alexandria, VA 22333

HQDA  
ATTN: DASG-PSP-E  
Falls Church, VA 22041-3258

HQDA  
ATTN: DAEN-RDM  
20 Massachusetts, NW  
Washington, D.C. 20314

END

DATE

FILMED

DTIC

JULY 88